LavaLAMP[™] DNA Master Mix

Lucigen[®] Simplifying Genomics

(CAT. NO. 30066-1, 30067-1)

- 1. Thaw all kit components on ice and set up reactions on ice.
- 2. Mix each component thoroughly before use by vortexing for 3 10 sec. Spin down contents at 4°C.
- **3.** Prepare reaction mixes in the order listed below in Table 1. When possible, make enough of the reaction mixes for all your reactions at one time.

Note: The reaction conditions recommended in Table 1 assume you will use the Green Fluorescent Dye available with Cat. No. 30067-1 and 1 μ L of Target DNA Sample in step 6. If you are not using the Green Fluorescent Dye and/or are using a different volume of Target DNA Sample, adjust the amount of nuclease-free H₂O accordingly.

	Table 1. F	Positive	Control, No	Target	Control	and	Experimental	Reaction	Setup
--	------------	----------	-------------	--------	---------	-----	--------------	----------	-------

	Positive Control	No Target Control (NTC)	Experimental
Component	Amount (μL)	Amount (μL)	Amount (μL)
Nuclease-free H ₂ O	8	8	8
LavaLAMP [™] DNA Master Mix (2X)	12.5	12.5	12.5
Target-Specific Primer Mix, 10X		2.5	2.5
DNA Positive Control LAMP Primer Mix	2.5		
Green Fluorescent Dye (optional)	1	1	1
Total volume	24	24	24

- 4. Mix the added reagents completely by pipetting.
- Dispense 24 μL each reaction mix into PCR tubes or a 96-well PCR plate.
 Note: In order to minimize cross-contamination, perform steps 6-8 in an area separate from area where you are preparing reaction mix.
- 6. Add 1 μL of Target DNA or Positive Control DNA to the appropriate reaction tubes or wells and 1 μL of nuclease-free water to the NTC reaction tubes or wells. Mix completely by pipetting.
- 7. Cap tubes or seal plate wells. Centrifuge briefly to collect contents prior to incubation.
- 8. Incubate the reactions as follows:

Step	Temperature	Time	
Optional: Preheat	90°C	≤5 min	
1. Amplification	Experimental and NTC : 68°C – 74°C Positive Control : 74°C	30 - 60 minutes	
2. Hold (Optional)	4°C	~	

- 9. Immediately stop amplification reactions using one of the three methods below. This step is required to stop the enzyme activity.
 - a. Hold on ice or at 4°C.
 - b. Add gel loading dye that produces a final 10 mM EDTA concentration.
 - c. Perform a heat-kill step in a thermocycler or heat block at 95°C for 5 minutes
- 10. Detect amplified product using your detection method of choice.

Lucigen Corporation Technical Support: (888) 575-9695 | (608) 831-9011 | <u>techsupport@lucigen.com</u>