



AMRESCO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

EZ Vision Two, DNA Dye & Buffer 6X

Code: N650

Chemical Formula:	Not Applicable	Manufacture Date:	(batch specific)
Molecular Weight:		Expiration/Reassay Date:	(batch specific)
CAS #:	Not Applicable	Storage:	
Appearance:		Grade:	BIOTECHNOLOGY GRADE
Dark green viscous liquid.			

Additional Information

TEST	SPECIFICATION	DISPOSITION
Abs.@595-615nm (1:500mL Water)	0.5 - 0.9	PASS
DNase	NONE	PASS
Electrophoresis (P/F)	PASS	PASS

Spec Set: N650BIO

Analysis may have been rounded to significant digits in specification limits.
Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____
Title:

Date Printed: 08/15/2016
Page 1 of 1