



ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

Corporate Headquarters

AMRESKO, LLC

28600 Fountain Parkway
Solon, OH 44139 USA
TEL: 440-349-1199 FAX: 440-349-1182
WWW.AMRESKO-INC.COM
Email: info@amresco-inc.com

CERTIFICATE OF ANALYSIS

EZ-VISION BLUELIGHT DNA DYE

Code: 1B1680
Batch Number: 1906C190

Chemical Formula: Not Applicable Manufacture Date: 06/17/2016
Molecular Weight: Not applicable Reassay Date: 06/30/2017
CAS #: Not Applicable
Appearance: Storage: Room Temperature
Clear, dark orange solution. Grade: BIOTECHNOLOGY GRADE

TEST	SPECIFICATION	ANALYSIS	DISPOSITION
Abs. @Lambda Max	0.23 - 0.36	0.27	PASS
DNA Staining	Pass / Fail	Pass	PASS
DNase	NONE / Fail	NONE	PASS
Lambda Max	490 - 510 nm	499 nm	PASS

Spec Set: 1B1680BIO

Analysis may have been rounded to significant digits in specification limits.
Product meets analytical specifications of the grades listed.

Internal ID #: 546

Signature: 

Title: QC Manager

Date Printed: 07/06/2016

Page 1 of 1