



ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

AMRESKO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

Sodium Chloride

Code: 0241

Chemical Formula:	NaCl	Manufacture Date:	(batch specific)
Molecular Weight:	58.44	Expiration/Reassay Date:	(batch specific)
CAS #:	7647-14-5		
Appearance:		Storage:	Room Temperature
Fine white granular powder		Grade:	BIOTECHNOLOGY GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Calcium	$\leq 0.005 \%$	PASS
DNase	NONE	PASS
Heavy Metals	$\leq 0.0005 \%$	PASS
Iodide	$\leq 0.002 \%$	PASS
Iron	$\leq 0.0002 \%$	PASS
Loss on Drying	$\leq <0.5 \%$	PASS
Magnesium	$\leq 0.005 \%$	PASS
Nitrogen Compounds	$\leq 0.001 \%$	PASS
pH (5%, Water) @25C	5.0 - 9.0	PASS
Phosphate	$\leq 0.0005 \%$	PASS
Potassium	$\leq 0.005 \%$	PASS
Protease	NONE	PASS
Purity	$\geq 99.9 \%$	PASS
RNase	NONE	PASS
Sulfate	$\leq 0.0004 \%$	PASS

Spec Set: 0241BIO

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____

Title:

Date Printed:

02/29/2016

Page 1 of 1