



AMRESKO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

RNase Inhibitor

Code: E633

Chemical Formula:	Not Applicable	Manufacture Date:	(batch specific)
Molecular Weight:		Expiration/Reassay Date:	(batch specific)
CAS #:	Not Applicable		
Appearance:		Storage:	
Light yellow liquid		Grade:	BIOTECHNOLOGY GRADE
		Additional Product Description:	
			from human placenta Solution in 20mM HEPES-KOH (pH 7.6), 50mM KCl, 5mM DDT, 50% Glycerol

Additional Information

TEST	SPECIFICATION	DISPOSITION
Activity @25C	>= 10 U/ul	PASS
DNase (endonuclease)	NONE	PASS
DNase (exonuclease)	NONE	PASS

Spec Set: E633BIO

Analysis may have been rounded to significant digits in specification limits.
Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____

Title:

Date Printed: 06/15/2016

Page 1 of 1