



AMRESCO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

DMSO (Dimethyl Sulfoxide)

Code: 0231

Chemical Formula:	(CH ₃) ₂ SO	Manufacture Date:	(batch specific)
Molecular Weight:	78.13	Expiration/Reassay Date:	(batch specific)
CAS #:	67-68-5		
Appearance:		Storage:	Room Temperature
Clear, colorless liquid		Grade:	ACS GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Abs.@270nm	<= 0.40	PASS
Abs.@290nm	<= 0.18	PASS
Abs.@310nm	<= 0.06	PASS
Abs.@330nm	<= 0.02	PASS
Abs.@350-400nm	<= 0.01	PASS
Purity	>= 99.9 %	PASS
Residue After Evaporation	<= 0.01 %	PASS
Titration Free Acid	<= 0.001 meq/g	PASS
Water (KF)	<= 0.1 %	PASS

Spec Set: 0231ACS

Analysis may have been rounded to significant digits in specification limits.
Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____

Title:

Date Printed: 02/29/2016