

Distributor

AMRESCO INC.

6681 Cochran Road

SOLON, OHIO 44139

VWR
(amresco) LIFE SCIENCE
Issuing Date 01-May-2009

Revision Date 08-Jan-2016

Revision Number 3

1. PRODUCT AND COMPANY IDENTIFICATION

Product Identifier Product Name

Dimethyl Sulfoxide (DMSO)

Other means of identification Product Code(s) Synonyms

0231 No information available

Recommended use of the chemical and restrictions on useRecommended UseFor Further Manufacturing Use Only.Uses advised againstNot for Human or Animal Drug Use

Details of the supplier of the safety data sheet

Company Address VWR International, LLC Radnor Corporate Center 100 Matsonford Road Radnor, PA 19087-8660 Company Phone Number E-mail Address Manufacturer Address VWR International, LLC 28600 Fountain Parkway Solon, Ohio 44139

1.800.829.2805 info@amresco-inc.com

Chemtrec 1-800-424-9300

Emergency Telephone Number Emergency Telephone Number

2. HAZARDS IDENTIFICATION

Classification

OSHA Regulatory Status

This chemical is considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200)

Skin Corrosion/Irritation	Category 2
Serious Eye Damage/Eye Irritation	Category 2
Specific target organ toxicity (single exposure)	Category 3
Flammable liquids	Category 4

Label elements

Emergency Overview

Warning

Hazard statements

May cause respiratory irritation. May cause drowsiness or dizziness H315 - Causes skin irritation H319 - Causes serious eye irritation H227 - Combustible liquid



I

Appearance Clear Colorless	Physical State Liquid	Odor No information available
----------------------------	-----------------------	-------------------------------

Precautionary Statements - Prevention

L

Wash face, hands and any exposed skin thoroughly after handling Avoid breathing dust/fume/gas/mist/vapors/spray Use only outdoors or in a well-ventilated area Keep cool

Precautionary Statements - Response

Specific treatment (see .? on this label) IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing If eye irritation persists: Get medical advice/attention IF ON SKIN: Wash with plenty of water and soap If skin irritation occurs: Get medical advice/attention Take off contaminated clothing and wash it before reuse IF INHALED: Remove person to fresh air and keep comfortable for breathing Call a POISON CENTER or doctor if you feel unwell In case of fire: Use CO2, dry chemical, or foam to extinguish

Precautionary Statements - Storage

Store in a well-ventilated place. Keep container tightly closed Store locked up

Precautionary Statements - Disposal

Dispose of contents/container to an approved waste disposal plant

Hazards not otherwise classified (HNOC) Not regulated Other Information

3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical Name	CAS-No	EC No.	Weight %	Trade Secret
Dimethyl Sulfoxide	67-68-5	200-664-3	95-100	Not applicable

4. FIRST AID MEASURES

First Aid Measures

Eye Contact	Rinse thoroughly with plenty of water for at least 15 minutes, lifting lower and upper eyelids. Consult a physician.	
Skin Contact	Readily permeates skin: substances mixed with DMSO may be quickly absorbed. Wash off immediately with soap and plenty of water removing all contaminated clothes and shoes.	
Inhalation	Move to fresh air in case of accidental inhalation of vapors or decomposition products. Apply artificial respiration if victim is not breathing.	
Ingestion	Never give anything by mouth to an unconscious person. Rinse mouth.	
Most important symptoms and effects, both acute and delayed		
Most Important Symptoms/Effects	No information available.	
Indication of any immediate medical attention and special treatment needed		

Notes to Physician Treat symptomatically.

5. FIRE-FIGHTING MEASURES

Suitable Extinguishing Media

Dry chemical, CO2, water spray or regular foam.

Unsuitable Extinguishing Media No information available.

Specific Hazards Arising from the Chemical

No information available.

Hazardous Combustion Carbon oxides. Sulfur oxides. Products

Explosion Data Sensitivity to Mechanical Impact None. Sensitivity to Static Discharge None.

Protective Equipment and Precautions for Firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

	6. ACCIDENTAL RELEASE MEASURES	
Personal precautions, protective equipment and emergency procedures		
Personal protection	Ensure adequate ventilation, especially in confined areas.	
Environmental Precautions		
Environmental Precautions	See Section 12 for additional Ecological information.	
Methods and material for containment and cleaning up		
Methods for Containment	Prevent further leakage or spillage if safe to do so.	
Methods for Cleaning Up	Absorb on sand or vermiculite and place in closed containers for disposal. Ventilate area and wash spill site after material pickup is complete. Use non-sparking tools and equipment.	
7. HANDLING AND STORAGE		

Precautions for Safe Handling

Handling	Handle in accordance with good industrial hygiene and safety practice. Hygroscopic. Light Sensitive. Heat Sensitive.
Conditions for safe storage, includ	ing any incompatibilities
Storage	Keep containers tightly closed in a dry, cool and well-ventilated place. Hygroscopic. Protect from light and moisture.
Incompatible Products	Acid chlorides. Strong acids. Strong oxidizing agents. Strong reducing agents. Phosphorus halides.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Control parameters

Exposure Guidelines

This product does not contain any hazardous materials with occupational exposure limits established by the region specific regulatory bodies.

Chemical Name	ACGIH TLV	OSHA PEL	NIOSH IDLH
Dimethyl Sulfoxide	-	-	-

67-68-5	
Appropriate engineering controls	
Engineering Measures	Showers Eyewash stations Ventilation systems.
Individual protection measures,	uch as personal protective equipment
Eye/Face Protection	Safety glasses with side-shields.
Skin and Body Protection	Wear protective gloves/clothing.

If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be provided in accordance with current local regulations.

Hygiene Measures Handle in accordance with good industrial hygiene and safety practice.

9. PHYSICAL AND CHEMICAL PROPERTIES

Information on basic physical and chemical properties

Physical State Appearance Color	Liquid Clear Colorless No information available	Odor Odor Threshold	No information available No information available
Property pH Melting point/freezing point Boiling Point/Range Flash Point (High in °C) Evaporation Rate Flammability (solid, gas) Flammability Limit in Air Upper flammability limit: Lower flammability limit: Lower flammability limit: Vapor pressure Vapor Density Specific Gravity Water Solubility Solubility in other solvents Partition coefficient Autoignition Temperature Decomposition Temperature Kinematic viscosity Dynamic viscosity Explosive Properties Oxidizing Properties	Values No information available 18 °C 189 °C 89 °C / 192 °F No information available No information	<u>Remarks • Method</u>	
Other Information Softening Point Molecular Weight	No information available No information available		
VOC Content Density Bulk Density	No information available 2.7 No information available		

10. STABILITY AND REACTIVITY

Reactivity

No data available

Chemical Stability Stable under recommended storage conditions. Possibility of Hazardous Reactions None under normal processing. Conditions to Avoid Hygroscopic. Protect from heat and light. Incompatible Materials Acid chlorides. Strong acids. Strong oxidizing agents. Strong reducing agents. Phosphorus halides. Hazardous Decomposition Products Carbon oxides. Sulfur oxides.

11. TOXICOLOGICAL INFORMATION

Information on likely routes of exposure

Product Information

Inhalation	May cause drownsiness and dizziness based on components. May cause irritation of respiratory tract.
Eye Contact	Irritating to eyes.
Skin Contact	Irritating to skin.
Ingestion	There is no data available for this product.

Chemical Name	LD50 Oral	LD50 Dermal	LC50 Inhalation
Dimethyl Sulfoxide	= 14500 mg/kg (Rat)	= 40 g/kg (Rat)	-
67-68-5			

Information on toxicological effects

Symptoms

No information available.

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Sensitization Mutagenic Effects Carcinogenicity	No informatio No informatio No informatio	on available.		
Chemical Name	ACGIH	IARC	NTP	OSHA
Dimethyl Sulfoxide 67-68-5	-	-	-	-
Reproductive Toxicity STOT - single exposure STOT - repeated exposur Target Organ Effects Aspiration hazard	No information No information Eyes, Skin.	No information available. No information available. No information available. Eyes, Skin. No information available.		

Numerical measures of toxicity - Product Information

The following values are calculated based on chapter 3.1 of the GHS document .

12. ECOLOGICAL INFORMATION

Ecotoxicity

100% of the mixture consists of components(s) of unknown hazards to the aquatic environment

Chemical Name	Algae/aquatic plants	Toxicity to Fish	Crustacea
Dimethyl Sulfoxide	12350 - 25500: 96 h Skeletonema	33 - 37: 96 h Oncorhynchus mykiss	7000: 24 h Daphnia species mg/L
67-68-5	costatum mg/L EC50	g/L LC50 static 41.7: 96 h Cyprinus	ÉC50

carpio g/L LC50 34000: 96 l Pimephales promelas mg/L LC 40: 96 h Lepomis macrochirus	50
LC50 static	g/ L

Persistence and Degradability

No information available.

Bioaccumulation/Accumulation

No information available.

Chemical Name	Partition coefficient
Dimethyl Sulfoxide	-
67-68-5	

Other Adverse Effects

No information available

13. DISPOSAL CONSIDERATIONS

Waste treatment methods

Waste Disposal Method Dispose of material in accordance with all federal, state, and local regulations.

Contaminated Packaging

Do not re-use empty containers.

Chemical Name	RCRA - Halogenated Organic Compounds	RCRA - P Series Wastes	RCRA - F Series Wastes	RCRA - K Series Wastes
Dimethyl Sulfoxide 67-68-5	-	-	-	-

Chemical Name	California Hazardous Waste Status
Dimethyl Sulfoxide	-
67-68-5	

14. TRANSPORT INFORMATION

DOT

Not regulated

IATA Not regulated

15. REGULATORY INFORMATION

International Inventories	
TSCA	Complies
DSL/NDSL	Complies
EINECS/ELINCS	Complies
ENCS	Complies
IECSC	Complies
KECL	Complies
PICCS	Complies
AICS	Complies

Legend:

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS - European Inventory of Existing Commercial Chemical Substances/EU List of Notified Chemical Substances

ENCS - Japan Existing and New Chemical Substances

IECSC - China Inventory of Existing Chemical Substances

KECL - Korean Existing and Evaluated Chemical Substances

PICCS - Philippines Inventory of Chemicals and Chemical Substances

AICS - Australian Inventory of Chemical Substances

U.S. Federal Regulations

SARA 313

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

Component	SARA 313 - Threshold Values %
Dimethyl Sulfoxide	-
67-68-5 (95-100)	

SARA 311/312 Hazard Categories

Acute Health Hazard	Yes
Chronic Health Hazard	No
Fire Hazard	Yes
Sudden Release of Pressure Hazard	No
Reactive Hazard	No

Clean Water Act

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42)

Chemical Name	CWA - Reportable Quantities	CWA - Toxic Pollutants	CWA - Priority Pollutants	CWA - Hazardous Substances
Dimethyl Sulfoxide 67-68-5	-	-	-	-

CERCLA

This material, as supplied, does not contain any substances regulated as hazardous substances under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302) or the Superfund Amendments and Reauthorization Act (SARA) (40 CFR 355). There may be specific reporting requirements at the local, regional, or state level pertaining to releases of this material

Chemical Name	Hazardous Substances RQs	Extremely Hazardous Substances RQs	RQ
Dimethyl Sulfoxide 67-68-5	-	-	-

U.S. State Regulations

California Proposition 65

This product does not contain any Proposition 65 chemicals.

Chemical Name	California Prop. 65
Dimethyl Sulfoxide - 67-68-5	-
U.S. State Pight to Know Pegulations	

U.S. State Right-to-Know Regulations

Chemical Name	New Jersey	Massachusetts	Pennsylvania
Dimethyl Sulfoxide	Х	-	-
67-68-5			

U.S. EPA Label Information

EPA Pesticide Registration Number Not regulated

16. OTHER INFORMATION

01-May-2009 08-Jan-2016

Issuing Date
Revision Date
Revision Note
No information available
Disclaimer
The chave information is half

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with

regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. VWR International and its Affiliates shall not be held liable for any damage resulting from handling.

End of Safety Data Sheet