



ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

AMRESKO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

Isopropyl Alcohol

Code: 0918

| | | | |
|------------------------|----------------------------------|--------------------------|---------------------|
| Chemical Formula: | C ₃ H ₇ OH | Manufacture Date: | (batch specific) |
| Molecular Weight: | 60.10 | Expiration/Reassay Date: | (batch specific) |
| CAS #: | 67-63-0 | | |
| Appearance: | | Storage: | Room Temperature |
| Clear colorless liquid | | Grade: | BIOTECHNOLOGY GRADE |

Additional Information

| TEST | SPECIFICATION | DISPOSITION |
|---------------|--------------------|-------------|
| Color (APHA) | <= 10 | PASS |
| Density | 0.782 - 0.788 g/mL | PASS |
| Moisture (KF) | <= 0.2 % | PASS |
| Purity | >= 99.0 % | PASS |

Spec Set: 0918BIO

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____

Title:

Date Printed: 06/09/2016

Page 1 of 1