



AMRESKO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

PBS 10X Solution

Code: J373

Chemical Formula:	N/A	Manufacture Date:	(batch specific)
Molecular Weight:		Expiration/Reassay Date:	(batch specific)
CAS #:	N/A		
Appearance:		Storage:	
Clear, colorless, haze-free liquid		Grade:	ULTRA PURE GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Conductivity (10%, Water) @25C	14000 - 17800 umhos	PASS
DNase	NONE	PASS
pH (10%, Water) @25C	7.3 - 7.5	PASS
pH @25C	6.7 - 6.9	PASS
Protease	NONE	PASS
RNase	NONE	PASS

Spec Set: J373UP

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____

Title:

Date Printed:

06/20/2016

Page 1 of 1