

Deoxynucleotide Triphosphates (dNTP) Pre-Mixed

Catalog number: PT-8014

Quantity: 1 ml

Concentration: 2.5 mM for each of dATP, dCTP, dGTP and dTTP

Description:

The dNTP is made by dissolving ultra-pure triphosphate powders in buffer at pH 7.4. The concentration was quantitated by spectral analysis and then adjusted to a final concentration of 2.5 mM.

Storage condition:

The solutions are stable at -20°C. However, for long term storage, we recommend storing in aliquots at -70°C.

*****Research Use Only*****

Please do not hesitate to contact us while you have any questions

Manufactured for and distributed by Protech Technology Enterprise Co.,Ltd

TEL : +886-2-2655-7677 / FAX : +886-2-2655-7601 / Toll Free : 0800-231-530

E-MAIL: service@bio-protech.com.tw ; tech@bio-protech.com.tw

WEB: <http://www.bio-protech.com.tw>