

SAFETY DATA SHEET

Issuing Date 19-May-2009 Revision Date 08-Jan-2016 Revision Number 2

1. PRODUCT AND COMPANY IDENTIFICATION

Product Identifier

Product Name Ethidium Bromide Stock Solution

Other means of identification

Product Code(s) X328 UN/ID no. 2810

Synonyms No information available

Recommended use of the chemical and restrictions on use

Recommended Use For Further Manufacturing Use Only.
Uses advised against Not for Human or Animal Drug Use

Details of the supplier of the safety data sheet

Company AddressManufacturer AddressVWR International, LLCVWR International, LLCRadnor Corporate Center28600 Fountain Parkway100 Matsonford RoadSolon, Ohio 44139

Radnor, PA 19087-8660

Company Phone Number 1-800-448-4442 E-mail Address info@amresco-inc.com

Emergency Telephone Number

Emergency Telephone Number Chemtrec 1-800-424-9300

2. HAZARDS IDENTIFICATION

Classification

OSHA Regulatory Status

This chemical is considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200)

Germ cell mutagenicity Category 2

Label elements

Emergency Overview

Warning

Hazard statements

H341 - Suspected of causing genetic defects



Appearance Dark red haze-free Physical State Liquid Odor No information available

Precautionary Statements - Prevention

Obtain special instructions before use

Do not handle until all safety precautions have been read and understood Use personal protective equipment as required

Precautionary Statements - Response

IF exposed or concerned: Get medical advice/attention

Precautionary Statements - Storage

Store locked up

Precautionary Statements - Disposal

Dispose of contents/container to an approved waste disposal plant

Hazards not otherwise classified (HNOC)

Not regulated

Other Information

Unknown Acute Toxicity 0 % of the mixture consists of ingredient(s) of unknown toxicity

3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical Name	CAS-No	EC No.	Weight %	Trade Secret
3,8-Diamino-1-ethyl-6-ph	1239-45-8	214-984-6	1	Not applicable
enylphenantridinium				
bromide				

4. FIRST AID MEASURES

First Aid Measures

Eye Contact Rinse thoroughly with plenty of water for at least 15 minutes, lifting lower and upper eyelids.

Consult a physician.

Skin ContactWash off immediately with soap and plenty of water removing all contaminated clothes and

shoes.

Inhalation Move to fresh air. If breathing becomes difficult, give oxygen. Consult a physician.

Ingestion Never give anything by mouth to an unconscious person. Rinse mouth.

Most important symptoms and effects, both acute and delayed

Most Important Symptoms/Effects No information available.

Indication of any immediate medical attention and special treatment needed

Notes to Physician Treat symptomatically.

5. FIRE-FIGHTING MEASURES

Suitable Extinguishing Media

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Unsuitable Extinguishing Media No information available.

Specific Hazards Arising from the Chemical

No information available.

Hazardous Combustion

Products

Emits toxic fumes.

Explosion Data

Sensitivity to Mechanical Impact None. Sensitivity to Static Discharge None.

Protective Equipment and Precautions for Firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures

Personal protection Ensure adequate ventilation, especially in confined areas.

Environmental Precautions

Environmental Precautions See Section 12 for additional Ecological information.

Methods and material for containment and cleaning up

Methods for Containment Prevent further leakage or spillage if safe to do so.

Methods for Cleaning Up Absorb with sand or vermiculite. Pick up and transfer to properly labeled containers.

Ventilate area and wash spill site after material pickup is complete.

7. HANDLING AND STORAGE

Precautions for Safe Handling

Handling Handle in accordance with good industrial hygiene and safety practice.

Conditions for safe storage, including any incompatibilities

Storage Keep containers tightly closed in a dry, cool and well-ventilated place.

Incompatible Products Strong oxidizing agents.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Control parameters

Exposure Guidelines This product does not contain any hazardous materials with occupational exposure limits

established by the region specific regulatory bodies.

Chemical Name	ACGIH TLV	OSHA PEL	NIOSH IDLH
3,8-Diamino-1-ethyl-6-phenylphena	-	-	-
ntridinium bromide			
1239-45-8			

Appropriate engineering controls

Engineering Measures Showers

Eyewash stations Ventilation systems.

Individual protection measures, such as personal protective equipment

Eye/Face Protection Safety glasses with side-shields. Avoid contact with eyes.

Skin and Body Protection Wear protective gloves/clothing.

Respiratory Protection If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved

respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be provided in accordance with current local regulations.

Remarks • Method

Hygiene Measures Handle in accordance with good industrial hygiene and safety practice.

9. PHYSICAL AND CHEMICAL PROPERTIES

Information on basic physical and chemical properties

Physical State Liquid

Appearance Dark red haze-free Odor No information available Color No information available **Odor Threshold** No information available

Property Values

No information available Hq

Melting point/freezing point

Boiling Point/Range No information available Flash Point (High in °C) No information available **Evaporation Rate** No information available Flammability (solid, gas) No information available

Flammability Limit in Air

Upper flammability limit: No information available Lower flammability limit: No information available Vapor pressure No information available Vapor Density No information available **Specific Gravity** No information available Water Solubility No information available Solubility in other solvents No information available No information available **Partition coefficient Autoignition Temperature** No information available **Decomposition Temperature** No information available Kinematic viscosity No information available Dynamic viscosity No information available **Explosive Properties** No information available **Oxidizing Properties** No information available

Other Information

No information available **Softening Point** Molecular Weight No information available **VOC Content** No information available **Density** No information available **Bulk Density** No information available

10. STABILITY AND REACTIVITY

Reactivity

No data available

Chemical Stability

Stable under recommended storage conditions.

Possibility of Hazardous Reactions

None under normal processing.

Conditions to Avoid

Extremes of temperature and direct sunlight.

Incompatible Materials

Strong oxidizing agents.

Hazardous Decomposition Products

Carbon oxides. Nitrogen oxides (NOx). Hydrogen bromide.

11. TOXICOLOGICAL INFORMATION

Information on likely routes of exposure

Product Information

Inhalation Harmful by inhalation.

Eye Contact There is no data available for this product.

Skin Contact There is no data available for this product.

Ingestion Harmful if swallowed.

Chemical Name	LD50 Oral	LD50 Dermal	LC50 Inhalation
3,8-Diamino-1-ethyl-6-phenylphena	=	-	-
ntridinium bromide			
1239-45-8			

Information on toxicological effects

Symptoms No information available.

Delayed and immediate effects as well as chronic effects from short and long-term exposure

SensitizationNo information available.Mutagenic EffectsNo information available.CarcinogenicityNo information available.

Chemical Name	ACGIH	IARC	NTP	OSHA
3,8-Diamino-1-ethyl-6-pheny	-	-	-	-
lphenantridinium bromide				
1230-45-8				

Reproductive Toxicity
STOT - single exposure
STOT - repeated exposure
No information available.
No information available.

Chronic Toxicity May cause heritable genetic damage.

Aspiration hazard No information available.

Numerical measures of toxicity - Product Information

Unknown Acute Toxicity 0 % of the mixture consists of ingredient(s) of unknown toxicity

The following values are calculated based on chapter 3.1 of the GHS document

12. ECOLOGICAL INFORMATION

Ecotoxicity

1% of the mixture consists of components(s) of unknown hazards to the aquatic environment

. 70 0. 10	periorite(e) er armiremi ma_arae	to the diquate of the original	
Chemical Name	Algae/aquatic plants	Toxicity to Fish	Crustacea
3,8-Diamino-1-ethyl-6-phenylphena	-	-	-
ntridinium bromide			
1239-45-8			

Persistence and Degradability

No information available.

Bioaccumulation/Accumulation

No information available.

Chemical Name	Partition coefficient
3,8-Diamino-1-ethyl-6-phenylphenantridinium bromide	-

1239-45-8

Other Adverse Effects No information available

13. DISPOSAL CONSIDERATIONS

Waste treatment methods

Waste Disposal Method Dispose of material in accordance with all federal, state, and local regulations.

Contaminated Packaging Dispose of in accordance with all federal, state and local regulations.

Chemical Name	RCRA	RCRA - Basis for Listing	RCRA - D Series Wastes	RCRA - U Series Wastes
3,8-Diamino-1-ethyl-6-pheny	-	-	-	-
Iphenantridinium bromide				
1239-45-8				

Chemical Name	RCRA - Halogenated Organic Compounds	RCRA - P Series Wastes	RCRA - F Series Wastes	RCRA - K Series Wastes
3,8-Diamino-1-ethyl-6-pheny	-	-	=	-
Iphenantridinium bromide				
1239-45-8				

Chemical Name	California Hazardous Waste Status
3,8-Diamino-1-ethyl-6-phenylphenantridinium bromide	-
1239-45-8	

14. TRANSPORT INFORMATION

DOT

UN/ID no. 2810

Proper shipping name Toxic, liquid, organic, N.O.S. (contains Ethidium Bromide)

Hazard Class 6.1 Packing Group III

IATA

UN/ID no. 2810

Proper shipping name Toxic, liquid, organic, N.O.S. (contains Ethidium Bromide)

Hazard Class 6.1 Packing Group III

15. REGULATORY INFORMATION

International Inventories

TSCA Does not comply
DSL/NDSL Does not comply
EINECS/ELINCS Complies

ENCS Does not Comply

IECSCCompliesKECLCompliesPICCSCompliesAICSComplies

Legend:

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS - European Inventory of Existing Commercial Chemical Substances/EU List of Notified Chemical Substances

ENCS - Japan Existing and New Chemical Substances IECSC - China Inventory of Existing Chemical Substances KECL - Korean Existing and Evaluated Chemical Substances

PICCS - Philippines Inventory of Chemicals and Chemical Substances

AICS - Australian Inventory of Chemical Substances

U.S. Federal Regulations

SARA 313

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

Component	SARA 313 - Threshold Values %
3,8-Diamino-1-ethyl-6-phenylphenantridinium bromide	-
1239-45-8 (1)	

SARA 311/312 Hazard Categories

Acute Health Hazard	Yes
Chronic Health Hazard	Yes
Fire Hazard	No
Sudden Release of Pressure Hazard	No
Reactive Hazard	No

Clean Water Act

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42)

Chemical Name	CWA - Reportable Quantities	CWA - Toxic Pollutants	CWA - Priority Pollutants	CWA - Hazardous Substances
3,8-Diamino-1-ethyl-6-pheny	-	=	-	-
Iphenantridinium bromide				
1239-45-8				

CERCLA

This material, as supplied, does not contain any substances regulated as hazardous substances under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302) or the Superfund Amendments and Reauthorization Act (SARA) (40 CFR 355). There may be specific reporting requirements at the local, regional, or state level pertaining to releases of this material

Chemical Name	Hazardous Substances RQs	Extremely Hazardous Substances RQs	RQ
3,8-Diamino-1-ethyl-6-phenylphena	•	-	-
ntridinium bromide			
1239-45-8			

U.S. State Regulations

California Proposition 65

This product does not contain any Proposition 65 chemicals.

Chemical Name	California Prop. 65	
3,8-Diamino-1-ethyl-6-phenylphenantridinium bromide - 1239-45-8	-	

U.S. State Right-to-Know Regulations

Chemical Name	New Jersey	Massachusetts	Pennsylvania
3,8-Diamino-1-ethyl-6-phenylphena	-	-	-
ntridinium bromide			
1239-45-8			

U.S. EPA Label Information

EPA Pesticide Registration Number Not regulated

4.0	OTHER	NICODM	ATION
1 h	()IHERI	$M \vdash C) \bowtie M$	Δ I I Ω

Issuing Date 19-May-2009
Revision Date 08-Jan-2016
Revision Note

Revision Date 08-Jan-2016

No information available

Disclaimer

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. VWR International and its Affiliates shall not be held liable for any damage resulting from handling.

End of Safety Data Sheet