



**AMRESKO LLC**

28600 Fountain Parkway  
Solon, Ohio USA 44139  
440/349-1199 FAX: 440/349-1182  
www.amresco-inc.com  
Email: info@amresco-inc.com

ISO 9001 CERTIFIED  
ISO 13485 CERTIFIED

## CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

### RIPA Lysis Buffer

**Code: N653**

Chemical Formula:	Not Applicable	Manufacture Date:	(batch specific)
Molecular Weight:		Expiration/Reassay Date:	(batch specific)
CAS #:	n/a		
Appearance:		Storage:	
Clear, colorless solution.		Grade:	BIOTECHNOLOGY GRADE

### Additional Information

TEST	SPECIFICATION	DISPOSITION
pH @25C	7.4 - 7.6	PASS
Tris Titration	48 - 55 mM - mMol	PASS

Spec Set: N653BIO

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: \_\_\_\_\_

Title:

Date Printed:

08/15/2016

Page 1 of 1