

# SAFETY DATA SHEET

Issuing Date 23-Jun-2010 Revision Date 08-Jan-2016 Revision Number 2

# 1. PRODUCT AND COMPANY IDENTIFICATION

**Product Identifier** 

Product Name Sodium Hydroxide (Pellets)

Other means of identification

Product Code(s) 0583 UN/ID no. 1823

**Synonyms** Caustic soda ; Lye ; Sodium hydrate ;

Recommended use of the chemical and restrictions on use

Recommended Use For Further Manufacturing Use Only.
Uses advised against Not for Human or Animal Drug Use

Details of the supplier of the safety data sheet

Company AddressManufacturer AddressVWR International, LLCVWR International, LLCRadnor Corporate Center28600 Fountain Parkway100 Matsonford RoadSolon, Ohio 44139

Radnor, PA 19087-8660

Company Phone Number 1-800-448-4442 E-mail Address info@amresco-inc.com

**Emergency Telephone Number** 

Emergency Telephone Number Chemtrec 1-800-424-9300

# 2. HAZARDS IDENTIFICATION

### Classification

# **OSHA Regulatory Status**

This chemical is considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200)

Skin Corrosion/Irritation	Category 1 Sub-category
Serious Eye Damage/Eye Irritation	Category 1
Corrosive to metals	Category 1

### Label elements

# **Emergency Overview**

### Danger

#### Hazard statements

H314 - Causes severe skin burns and eye damage

H290 - May be corrosive to metals



Appearance No information available Physical State Pellets Odor No information available

\_\_\_\_\_

#### **Precautionary Statements - Prevention**

Do not breathe dust/fume/gas/mist/vapors/spray

Wash face, hands and any exposed skin thoroughly after handling

Wear protective gloves/protective clothing/eye protection/face protection

Keep only in original container

#### **Precautionary Statements - Response**

Immediately call a POISON CENTER or doctor/physician

Specific treatment (see .? on this label)

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing Immediately call a POISON CENTER or doctor/physician

IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower

Wash contaminated clothing before reuse

IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing

Immediately call a POISON CENTER or doctor/physician IF SWALLOWED: Rinse mouth. DO NOT induce vomiting

Absorb spillage to prevent material damage

### **Precautionary Statements - Storage**

Store locked up

Store in corrosive resistant/ .? container with a resistant inner liner

#### **Precautionary Statements - Disposal**

Dispose of contents/container to an approved waste disposal plant

#### Hazards not otherwise classified (HNOC)

Not regulated

# Other Information

- Harmful to aquatic life with long lasting effects
- · Harmful to aquatic life

Unknown Acute Toxicity 0 % of the mixture consists of ingredient(s) of unknown toxicity

# 3. COMPOSITION/INFORMATION ON INGREDIENTS

**Synonyms** Caustic soda ;, Lye ;, Sodium hydrate ;.

Chemical Name	CAS-No	EC No.	Weight %	Trade Secret
Sodium hydroxide	1310-73-2	215-185-5	95-100	Not applicable

# 4. FIRST AID MEASURES

#### **First Aid Measures**

General Advice Consult a physician. Move out of dangerous area. Show this safety data sheet to the doctor

in attendance.

**Eye Contact** Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.

Continue rinsing eyes during transport to hospital.

Skin Contact Call a physician immediately. Wash off immediately with soap and plenty of water removing

all contaminated clothes and shoes. In case of contact with substance, keep exposed skin areas immersed in water or covered with wet bandages until medical attention is received.

**Inhalation** Move to fresh air. If not breathing, give artificial respiration. Consult a physician.

**Ingestion** Do not induce vomiting. Never give anything by mouth to an unconscious person. Rinse

mouth. Consult a physician.

### Most important symptoms and effects, both acute and delayed

**Most Important Symptoms/Effects** The most important known symptoms and effects are described in sections 2 & 11.

Indication of any immediate medical attention and special treatment needed

Notes to Physician Treat symptomatically.

# 5. FIRE-FIGHTING MEASURES

### **Suitable Extinguishing Media**

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Unsuitable Extinguishing Media No information available.

# **Specific Hazards Arising from the Chemical**

No information available.

**Hazardous Combustion** 

Sodium oxides.

**Products** 

**Explosion Data** 

Sensitivity to Mechanical Impact None. Sensitivity to Static Discharge None.

### **Protective Equipment and Precautions for Firefighters**

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

# 6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures

Personal protection Ensure adequate ventilation, especially in confined areas.

Environmental Precautions

**Environmental Precautions** See Section 12 for additional Ecological information.

Methods and material for containment and cleaning up

**Methods for Containment** Prevent further leakage or spillage if safe to do so.

Methods for Cleaning Up Avoid dust formation. Pick up and transfer to properly labeled containers. Ventilate area

and wash spill site after material pickup is complete.

# 7. HANDLING AND STORAGE

**Precautions for Safe Handling** 

Handling Handle in accordance with good industrial hygiene and safety practice. Protect from

moisture.

Conditions for safe storage, including any incompatibilities

**Storage** Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from

moisture.

Incompatible Products Strong oxidizing agents. Acids. Organic materials. Chlorinated solvents. Absorbs CO2 from

air. Aluminium. Phosphorus. Tin /Tin oxides. Zinc.

### 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

#### Control parameters

Exposure Guidelines

Chemical Name	ACGIH TLV	OSHA PEL	NIOSH IDLH
Sodium hydroxide	-	TWA: 2 mg/m <sup>3</sup>	IDLH: 10 mg/m <sup>3</sup>
1310-73-2		_	Ceiling: 2 mg/m <sup>3</sup>

# **Appropriate engineering controls**

**Engineering Measures** Showers

Eyewash stations Ventilation systems.

### Individual protection measures, such as personal protective equipment

**Eye/Face Protection** Avoid contact with eyes. Safety glasses with side-shields.

**Skin and Body Protection** Wear protective gloves/clothing.

**Respiratory Protection** If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved

respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be

provided in accordance with current local regulations.

**Hygiene Measures**Handle in accordance with good industrial hygiene and safety practice.

### 9. PHYSICAL AND CHEMICAL PROPERTIES

#### Information on basic physical and chemical properties

Physical State Pellets

AppearanceNo information availableOdorNo information availableColorNo information availableOdor ThresholdNo information available

<u>Property</u> <u>Values</u> <u>Remarks • Method</u>

pH 13 - 14 Melting point/freezing point 318 °C

Boiling Point/Range 1390 °C / 2534 °F
Flash Point (High in °C) No information available
Evaporation Rate No information available
Flammability (solid, gas) No information available
Flammability Limit in Air

Upper flammability limit:
Lower flammability limit:
Vapor pressure
Vapor Density

No information available
No information available
No information available
No information available

Specific Gravity 2.13

Water Solubility No information available Solubility in other solvents No information available Partition coefficient No information available **Autoignition Temperature** No information available **Decomposition Temperature** No information available Kinematic viscosity No information available **Dvnamic viscosity** No information available **Explosive Properties** No information available **Oxidizing Properties** No information available

#### Other Information

Softening PointNo information availableMolecular WeightNo information availableVOC ContentNo information availableDensityNo information availableBulk DensityNo information available

# 10. STABILITY AND REACTIVITY

#### Reactivity

No data available

#### **Chemical Stability**

Stable under recommended storage conditions.

### **Possibility of Hazardous Reactions**

None under normal processing.

### **Conditions to Avoid**

Do not allow water to enter container,

Heat of solution is very high, and with limited amounts of water , violent boiling may occur. Protect from moisture.

#### **Incompatible Materials**

Strong oxidizing agents. Acids. Organic materials. Chlorinated solvents. Absorbs CO2 from air. Aluminium. Phosphorus. Tin /Tin oxides. Zinc.

### **Hazardous Decomposition Products**

In the event of fire: see section 5.

# 11. TOXICOLOGICAL INFORMATION

# Information on likely routes of exposure

#### **Product Information**

**Inhalation**Contact with moist mucous membranes of the respiratory system can cause caustic

condition resulting in burns.

**Eye Contact** Corrosive to the eyes and may cause severe damage including blindness.

**Skin Contact** Contact causes severe skin irritation and possible burns.

**Ingestion** Causes burns.

Chemical Name	LD50 Oral	LD50 Dermal	LC50 Inhalation
Sodium hydroxide	-	= 1350 mg/kg (Rabbit)	-
1310-73-2			

# Information on toxicological effects

**Symptoms** No information available.

# Delayed and immediate effects as well as chronic effects from short and long-term exposure

SensitizationNo information available.Mutagenic EffectsNo information available.CarcinogenicityNo information available.

Chemical Name	ACGIH	IARC	NTP	OSHA
Sodium hydroxide	-	-	-	-
1310-73-2				

Reproductive Toxicity
STOT - single exposure
STOT - repeated exposure
Aspiration hazard
No information available.
No information available.
No information available.

# Numerical measures of toxicity - Product Information

Unknown Acute Toxicity 0 % of the mixture consists of ingredient(s) of unknown toxicity

# 12. ECOLOGICAL INFORMATION

### **Ecotoxicity**

0% of the mixture consists of components(s) of unknown hazards to the aquatic environment

Chemical Name	Algae/aquatic plants	Toxicity to Fish	Crustacea
Sodium hydroxide	-	45.4: 96 h Oncorhynchus mykiss	-
1310-73-2		mg/L LC50 static	

# Persistence and Degradability

No information available.

### **Bioaccumulation/Accumulation**

No information available.

Chemical Name	Partition coefficient
Sodium hydroxide	-
1310-73-2	

Other Adverse Effects No information available

# 13. DISPOSAL CONSIDERATIONS

### Waste treatment methods

Waste Disposal Method Dispose of material in accordance with all federal, state, and local regulations.

**Contaminated Packaging** Do not re-use empty containers.

Chemical Name	RCRA - Halogenated Organic Compounds	RCRA - P Series Wastes	RCRA - F Series Wastes	RCRA - K Series Wastes
Sodium hydroxide 1310-73-2	-	-	-	-

Chemical Name	California Hazardous Waste Status
Sodium hydroxide	-
1310-73-2	

# 14. TRANSPORT INFORMATION

DOT

UN/ID no. 1823

Proper shipping name SODIUM HYDROXIDE, SOLID

Hazard Class 8
Packing Group

**IATA** 

UN/ID no. 1823

Proper shipping name SODIUM HYDROXIDE, SOLID

Hazard Class 8
Packing Group ||

# 15. REGULATORY INFORMATION

**International Inventories** 

TSCA Complies
DSL/NDSL Complies
EINECS/ELINCS Complies
ENCS Complies
IECSC Complies
KECL Complies

PICCS Complies AICS Complies

Legend:

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS - European Inventory of Existing Commercial Chemical Substances/EU List of Notified Chemical Substances

**ENCS** - Japan Existing and New Chemical Substances

IECSC - China Inventory of Existing Chemical Substances

**KECL** - Korean Existing and Evaluated Chemical Substances

PICCS - Philippines Inventory of Chemicals and Chemical Substances

AICS - Australian Inventory of Chemical Substances

# U.S. Federal Regulations

#### **SARA 313**

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

Component	SARA 313 - Threshold Values %
Sodium hydroxide	-
1310-73-2 ( 95-100 )	

#### SARA 311/312 Hazard Categories

Acute Health Hazard Yes
Chronic Health Hazard No
Fire Hazard No
Sudden Release of Pressure Hazard No
Reactive Hazard No

#### **Clean Water Act**

This product contains the following substances which are regulated pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42)

Chemical Name	CWA - Reportable Quantities	CWA - Toxic Pollutants	CWA - Priority Pollutants	CWA - Hazardous Substances
Sodium hydroxide 1310-73-2	-	-	-	-

# CERCLA

This material, as supplied, contains one or more substances regulated as a hazardous substance under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302)

Chemical Name	Hazardous Substances RQs	Extremely Hazardous Substances RQs	RQ
Sodium hydroxide 1310-73-2	1000 lb	-	RQ 1000 lb final RQ RQ 454 kg final RQ

# U.S. State Regulations

# **California Proposition 65**

This product does not contain any Proposition 65 chemicals.

Chemical Name	California Prop. 65	
Sodium hydroxide - 1310-73-2	-	

# U.S. State Right-to-Know Regulations

Chemical Name	New Jersey	Massachusetts	Pennsylvania
Sodium hydroxide	X	X	X
1310-73-2			

# **U.S. EPA Label Information**

EPA Pesticide Registration Number Not regulated

16. OTHER INFORMATION

**Issuing Date** 23-Jun-2010 **Revision Date** 08-Jan-2016 **Revision Note** 

No information available

**Disclaimer** 

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. VWR International and its Affiliates shall not be held liable for any damage resulting from handling.

**End of Safety Data Sheet**