



**AMRESCO LLC**

28600 Fountain Parkway  
Solon, Ohio USA 44139  
440/349-1199 FAX: 440/349-1182  
www.amresco-inc.com  
Email: info@amresco-inc.com

ISO 9001 CERTIFIED  
ISO 13485 CERTIFIED

**CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS**

**Sodium Hydroxide (Pellets)**

**Code: 0583**

Chemical Formula:	NaOH	Manufacture Date:	(batch specific)
Molecular Weight:	40.00	Expiration/Reassay Date:	(batch specific)
CAS #:	1310-73-2		
Appearance:		Storage:	Room Temperature
Pellets		Grade:	ACS GRADE

**Additional Information**

TEST	SPECIFICATION	DISPOSITION
Calcium	<= 0.005 %	PASS
Chloride	<= 0.005 %	PASS
Heavy Metals	<= 0.002 %	PASS
Iron	<= 0.001 %	PASS
Magnesium	<= 0.002 %	PASS
Mercury	<= 0.1 ppm	PASS
Nickel	<= 0.001 %	PASS
Nitrogen Compounds	<= 0.001 %	PASS
Phosphate	<= 0.001 %	PASS
Potassium	<= 0.02 %	PASS
Purity	>= 97.0 %	PASS
Sodium Carbonate	<= 1.0 %	PASS
Sulfate	<= 0.003 %	PASS

Spec Set: 0583ACS

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: \_\_\_\_\_

Title:

Date Printed: 05/04/2016

Page 1 of 1