



ISO 9001 CERTIFIED
 ISO 13485 CERTIFIED

AMRESCO LLC
 28600 Fountain Parkway
 Solon, Ohio USA 44139
 440/349-1199 FAX: 440/349-1182
 www.amresco-inc.com
 Email: info@amresco-inc.com

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

EGTA

Code: 0732

Chemical Formula:	C14H24N2O10	Manufacture Date:	(batch specific)
Molecular Weight:	380.35	Expiration/Reassay Date:	(batch specific)
CAS #:	67-42-5		
Appearance:		Storage:	
Fine white crystalline powder		Grade:	ULTRA PURE GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
DNase	NONE	PASS
Heavy Metals (as Pb)	<= 0.001 %	PASS
Identification (IR)	PASS	PASS
Loss on Drying	<= 1.0 %	PASS
Melting Range w/Decomposition	>= 238 Degrees C	PASS
Protease	NONE	PASS
Purity	>= 97.0 %	PASS
RNase	NONE	PASS

Spec Set: 0732UP

Analysis may have been rounded to significant digits in specification limits.
 Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____
 Title:

Date Printed: 05/18/2016
 Page 1 of 1