



AMRESKO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

Proteinase K Solution, 20mg/mL

Code: E195

Chemical Formula:	Not Applicable	Manufacture Date:	(batch specific)
Molecular Weight:		Expiration/Reassay Date:	(batch specific)
CAS #:	Not Applicable		
Appearance:		Storage:	Cold
Slightly yellow, haze-free liquid.		Grade:	BIOTECHNOLOGY GRADE
		Additional Product Description:	
		20 mg/ml Proteinase K in a pH adjusted solution (pH 8.0) that is 10 mM Tris, 1 mM Calcium Chloride & 30% Glycerol.	

Additional Information

TEST	SPECIFICATION	DISPOSITION
Activity (37 degrees C)	>= 600 U/ml	PASS
RNase	NONE	PASS

Spec Set: E195BIO

Analysis may have been rounded to significant digits in specification limits.
Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____

Title:

Date Printed: 06/10/2016

Page 1 of 1