



**AMRESKO LLC**

28600 Fountain Parkway  
Solon, Ohio USA 44139  
440/349-1199 FAX: 440/349-1182  
www.amresco-inc.com  
Email: info@amresco-inc.com

ISO 9001 CERTIFIED  
ISO 13485 CERTIFIED

**CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS**

**NEXT GEL, 15%**

**Code: M258**

Chemical Formula:	Not Applicable	Manufacture Date:	(batch specific)
Molecular Weight:		Expiration/Reassay Date:	(batch specific)
CAS #:	Not Applicable	Storage:	
Appearance:		Grade:	PROTEOMICS GRADE
Clear, colorless liquid.			

**Additional Information**

TEST	SPECIFICATION	DISPOSITION
Conductivity (1:20, Water)	430 - 550 uS	PASS
Electrophoresis (P/F)	PASS	PASS
pH @25C	6.8 - 7.2	PASS
Polymerization Time	<= 30 minutes	PASS
Protease	NONE	PASS

Spec Set: M258PTM

Analysis may have been rounded to significant digits in specification limits.  
Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: \_\_\_\_\_

Title:

Date Printed: 08/11/2016

Page 1 of 1