



AMRESKO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

Terrific Broth Powder

Code: J869

Chemical Formula:	Not Applicable	Manufacture Date:	(batch specific)
Molecular Weight:		Expiration/Reassay Date:	(batch specific)
CAS #:	Not Applicable		
Appearance:		Storage:	
Powder		Grade:	BIOTECHNOLOGY GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Conductivity (10%, Water) @25C	1000 - 3000 umhos	PASS
pH (47.6g/1000ml Water) @25C	7.0 - 7.4	PASS
Solubility after autoclaving: 47.6g/1000mL water	Pass	PASS
Supports Growth of Bacteria	PASS	PASS

Spec Set: J869BIO

Analysis may have been rounded to significant digits in specification limits.
Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____
Title:

Date Printed: 07/27/2016
Page 1 of 1