

SAFETY DATA SHEET

Issuing Date 21-Apr-2009 Revision Date 08-Jan-2016 Revision Number 2

1. PRODUCT AND COMPANY IDENTIFICATION

Product Identifier

Product Name Silver Nitrate

Other means of identification

Product Code(s) 0377 UN/ID no. 1493

Synonyms Nitric acid, silver salt;

Recommended use of the chemical and restrictions on use

Recommended Use For Further Manufacturing Use Only.
Uses advised against Not for Human or Animal Drug Use

Details of the supplier of the safety data sheet

Company AddressManufacturer AddressVWR International, LLCVWR International, LLCRadnor Corporate Center28600 Fountain Parkway100 Matsonford RoadSolon, Ohio 44139

Radnor, PA 19087-8660

Company Phone Number 1-800-448-4442 E-mail Address info@amresco-inc.com

Emergency Telephone Number

Emergency Telephone Number Chemtrec 1-800-424-9300

2. HAZARDS IDENTIFICATION

Classification

OSHA Regulatory Status

This chemical is considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200)

| Acute Toxicity - Oral | Category 4 |
|-----------------------------------|---------------------------|
| Skin Corrosion/Irritation | Category 1 Sub-category B |
| Serious Eye Damage/Eye Irritation | Category 1 |
| Oxidizing solids | Category 2 |

Label elements

Emergency Overview

Danger

Hazard statements

H302 - Harmful if swallowed

H314 - Causes severe skin burns and eye damage

H272 - May intensify fire; oxidizer



Revision Date 08-Jan-2016

Appearance Fine White Physical State Crystalline Powder Odor No information available

Precautionary Statements - Prevention

Wash face, hands and any exposed skin thoroughly after handling

Do not eat, drink or smoke when using this product

Do not breathe dust/fume/gas/mist/vapors/spray

Wear protective gloves/protective clothing/eye protection/face protection

Keep away from heat/sparks/open flames/hot surfaces. - No smoking

Keep/Store away from clothing/ combustible materials

Take any precaution to avoid mixing with combustibles - .?

Precautionary Statements - Response

Immediately call a POISON CENTER or doctor/physician

Specific treatment (see .? on this label)

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing Immediately call a POISON CENTER or doctor/physician

IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower

Wash contaminated clothing before reuse

IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing

Immediately call a POISON CENTER or doctor/physician

IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell

Rinse mouth

Do NOT induce vomiting

In case of fire: Use CO2, dry chemical, or foam for extinction

Precautionary Statements - Storage

Store locked up

Precautionary Statements - Disposal

Dispose of contents/container to an approved waste disposal plant

Hazards not otherwise classified (HNOC)

Not regulated

Other Information

- Very toxic to aquatic life with long lasting effects
- · Very toxic to aquatic life

Unknown Acute Toxicity

0 % of the mixture consists of ingredient(s) of unknown toxicity

3. COMPOSITION/INFORMATION ON INGREDIENTS

Synonyms Nitric acid, silver salt;.

| Chemical Name | CAS-No | EC No. | Weight % | Trade Secret |
|----------------|-----------|-----------|----------|----------------|
| Silver nitrate | 7761-88-8 | 231-853-9 | 95-100 | Not applicable |

4. FIRST AID MEASURES

First Aid Measures

Eye Contact Rinse thoroughly with plenty of water for at least 15 minutes, lifting lower and upper eyelids.

Consult a physician.

Skin Contact Wash off immediately with soap and plenty of water removing all contaminated clothes and

shoes. Consult a physician.

Inhalation Move to fresh air. If not breathing, give artificial respiration. Consult a physician.

Ingestion Do NOT induce vomiting. Never give anything by mouth to an unconscious person. Rinse

mouth. Consult a physician.

Most important symptoms and effects, both acute and delayed

Most Important Symptoms/Effects No information available.

Indication of any immediate medical attention and special treatment needed

Notes to Physician Treat symptomatically.

5. FIRE-FIGHTING MEASURES

Suitable Extinguishing Media

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Unsuitable Extinguishing Media No information available.

Specific Hazards Arising from the Chemical

No information available.

Hazardous Combustion

Nitrogen oxides (NOx). Silver/Silver oxide. Container explosion may occur under fire

conditions.

Explosion Data

Products

Sensitivity to Mechanical Impact None. Sensitivity to Static Discharge None.

Protective Equipment and Precautions for Firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures

Personal protection Ensure adequate ventilation, especially in confined areas.

Environmental Precautions

Environmental Precautions See Section 12 for additional Ecological information.

Methods and material for containment and cleaning up

Methods for Containment Prevent further leakage or spillage if safe to do so.

Methods for Cleaning Up Avoid dust formation. Pick up and transfer to properly labeled containers. Ventilate area

and wash spill site after material pickup is complete.

7. HANDLING AND STORAGE

Precautions for Safe Handling

Handling Light Sensitive. Handle in accordance with good industrial hygiene and safety practice.

Conditions for safe storage, including any incompatibilities

Storage Keep container tightly closed in a dry and well-ventilated place. Protect from light.

Incompatible Products Strong reducing agents. Strong bases. Alcohols. Ammonia. Magnesium.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Control parameters

Exposure Guidelines

| Chemical Name | ACGIH TLV | OSHA PEL | NIOSH IDLH |
|----------------|-----------|----------|--------------|
| Silver nitrate | = | - | - |
| 7761-88-8 | | | |

Appropriate engineering controls

Engineering Measures Showers

Eyewash stations Ventilation systems.

Individual protection measures, such as personal protective equipment

Eye/Face Protection Avoid contact with eyes. Safety glasses with side-shields.

Skin and Body Protection Wear protective gloves/clothing.

Respiratory Protection If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved

respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be

provided in accordance with current local regulations.

Hygiene Measures Handle in accordance with good industrial hygiene and safety practice.

9. PHYSICAL AND CHEMICAL PROPERTIES

Information on basic physical and chemical properties

Physical State Crystalline Powder

AppearanceFine WhiteOdorNo information availableColorNo information availableOdor ThresholdNo information available

Property Values Remarks • Method

pH No information available

Melting point/freezing point 212 °C / 414 °F

Boiling Point/Range 440 °C / 824 °F

Flash Point (High in °C) No information available

Evaporation Rate No information available

Flammability (solid, gas) No information available

Flammability Limit in Air

Upper flammability limit: No information available No information available Lower flammability limit: Vapor pressure No information available Vapor Density No information available **Specific Gravity** No information available Water Solubility No information available Solubility in other solvents No information available Partition coefficient No information available **Autoignition Temperature** No information available **Decomposition Temperature** No information available Kinematic viscosity No information available **Dynamic viscosity** No information available

Explosive PropertiesNo information available
Oxidizing Properties
Oxidizer

Other Information

Softening Point
Molecular Weight
VOC Content
No information available
No information available
No information available

Density 4.35

Bulk Density No information available

10. STABILITY AND REACTIVITY

Reactivity

No data available

Chemical Stability

Stable under recommended storage conditions.

Possibility of Hazardous Reactions

None under normal processing.

Conditions to Avoid

May decompose on exposure to light. Protect from light.

Incompatible Materials

Strong reducing agents. Strong bases. Alcohols. Ammonia. Magnesium.

Hazardous Decomposition Products

In the event of fire: see section 5.

11. TOXICOLOGICAL INFORMATION

Information on likely routes of exposure

Product Information

InhalationContact with moist mucous membranes of the respiratory system can cause caustic

condition resulting in burns.

Eye Contact Corrosive to the eyes and may cause severe damage including blindness.

Skin Contact Contact causes severe skin irritation and possible burns.

Ingestion Harmful if swallowed.

| Chemical Name | LD50 Oral | LD50 Dermal | LC50 Inhalation |
|----------------|--------------------|-------------|-----------------|
| Silver nitrate | = 1173 mg/kg (Rat) | - | - |
| 7761-88-8 | | | |

Information on toxicological effects

Symptoms No information available.

Delayed and immediate effects as well as chronic effects from short and long-term exposure

SensitizationNo information available.Mutagenic EffectsNo information available.CarcinogenicityNo information available.

| Chemical Name | ACGIH | IARC | NTP | OSHA |
|----------------|-------|----------|-----|------|
| Silver nitrate | - | Group 2A | - | - |
| 7761-88-8 | | | | |

Reproductive Toxicity
STOT - single exposure
STOT - repeated exposure
Target Organ Effects

No information available.
No information available.
Liver.

Aspiration hazard No information available.

Numerical measures of toxicity - Product Information

Unknown Acute Toxicity 0 % of the mixture consists of ingredient(s) of unknown toxicity

The following values are calculated based on chapter 3.1 of the GHS document

12. ECOLOGICAL INFORMATION

Ecotoxicity

0% of the mixture consists of components(s) of unknown hazards to the aquatic environment

| Chemical Name | Algae/aquatic plants | Toxicity to Fish | Crustacea |
|----------------|------------------------------------|------------------------------------|---------------------------------|
| Silver nitrate | - | 0.00512 - 0.00787: 96 h Poecilia | 0.0006: 48 h Daphnia magna mg/L |
| 7761-88-8 | | reticulata mg/L LC50 semi-static | EC50 0.0008 - 0.0011: 48 h |
| | | 0.009 - 0.02: 96 h Lepomis | Daphnia magna mg/L EC50 Static |
| | | macrochirus mg/L LC50 | 0.0008 - 0.001: 48 h Daphnia |
| | | flow-through 0.0242 - 0.0484: 96 h | magna mg/L EC50 Flow through |
| | | Lepomis macrochirus mg/L LC50 | |
| | | semi-static 0.0064 - 0.0106: 96 h | |
| | | Pimephales promelas mg/L LC50 | |
| | | semi-static 0.0027: 96 h Cyprinus | |
| | | carpio mg/L LC50 semi-static | |
| | 0.0075: 96 h Oncorhynchus mykiss | | |
| | mg/L LC50 semi-static 0.00452 - | | |
| | | 0.00638: 96 h Pimephales promelas | |
| | | mg/L LC50 flow-through 0.00181 - | |
| | | 0.00214: 96 h Pimephales promelas | |
| | mg/L LC50 static 0.00839 - 0.1802: | | |
| | 96 h Oncorhynchus mykiss mg/L | | |
| | | LC50 static 0.05 - 0.07: 96 h | |
| | | Lepomis macrochirus mg/L LC50 | |
| | | static 0.001339 - 0.001637: 96 h | |
| | | Oncorhynchus mykiss mg/L LC50 | |
| | | flow-through 0.009: 96 h | |
| | | Pimephales promelas mg/L LC50 | |

Persistence and Degradability

No information available.

Bioaccumulation/Accumulation

No information available.

| Chemical Name | Partition coefficient |
|----------------|-----------------------|
| Silver nitrate | - |
| 7761-88-8 | |

Other Adverse Effects No information available

13. DISPOSAL CONSIDERATIONS

Waste treatment methods

Waste Disposal Method Dispose of material in accordance with all federal, state, and local regulations.

Contaminated Packaging Do not re-use empty containers.

| Chemical Name | RCRA - Halogenated Organic Compounds | RCRA - P Series Wastes | RCRA - F Series Wastes | RCRA - K Series Wastes |
|-----------------------------|---|------------------------|------------------------|------------------------|
| Silver nitrate 7761-88-8 | - | - | - | - |

| Chemical Name | California Hazardous Waste Status |
|----------------|-----------------------------------|
| Silver nitrate | Toxic |
| 7761-88-8 | |

14. TRANSPORT INFORMATION

DOT

UN/ID no. 1493

Proper shipping name SILVER NITRATE

Hazard Class 5.1 Packing Group II

Reportable Quantity (RQ) 1 LB / 0.454 KG

Emergency Response Guide 140

Number

IATA

UN/ID no. 1493

Proper shipping name SILVER NITRATE

Hazard Class 5.1
Packing Group II
ERG Code 140

15. REGULATORY INFORMATION

International Inventories

TSCA Complies **DSL/NDSL** Complies **EINECS/ELINCS** Complies Complies **ENCS** Complies **IECSC** Complies **KECL PICCS** Complies **AICS** Complies

Legend:

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS - European Inventory of Existing Commercial Chemical Substances/EU List of Notified Chemical Substances

ENCS - Japan Existing and New Chemical Substances

IECSC - China Inventory of Existing Chemical Substances

KECL - Korean Existing and Evaluated Chemical Substances

PICCS - Philippines Inventory of Chemicals and Chemical Substances

AICS - Australian Inventory of Chemical Substances

U.S. Federal Regulations

SARA 313

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

| Component | SARA 313 - Threshold Values % |
|----------------------|-------------------------------|
| Silver nitrate | - |
| 7761-88-8 (95-100) | |

SARA 311/312 Hazard Categories

Acute Health Hazard Yes
Chronic Health Hazard No
Fire Hazard Yes
Sudden Release of Pressure Hazard No
Reactive Hazard No

Clean Water Act

This product contains the following substances which are regulated pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42)

| Chemical Name | CWA - Reportable Quantities | CWA - Toxic Pollutants | CWA - Priority Pollutants | CWA - Hazardous Substances |
|----------------|--------------------------------|------------------------|---------------------------|-------------------------------|
| Silver nitrate | 1 lb | X | - | X |
| 7761-88-8 | | | | |

CERCLA

This material, as supplied, does not contain any substances regulated as hazardous substances under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302) or the Superfund Amendments and Reauthorization Act (SARA) (40 CFR 355). There may be specific reporting requirements at the local, regional, or state level

pertaining to releases of this material

| Che | emical Name | Hazardous Substances RQs | Extremely Hazardous Substances RQs | RQ |
|-----|-----------------------------|--------------------------|------------------------------------|--|
| _ | Silver nitrate 7761-88-8 | 1 lb | - | RQ 1 lb final RQ RQ 0.454 kg final RQ |

U.S. State Regulations

California Proposition 65

This product does not contain any Proposition 65 chemicals.

| Chemical Name | California Prop. 65 | |
|----------------------------|---------------------|--|
| Silver nitrate - 7761-88-8 | - | |

U.S. State Right-to-Know Regulations

| Chemical Name | New Jersey | Massachusetts | Pennsylvania |
|----------------|------------|---------------|--------------|
| Silver nitrate | X | X | X |
| 7761-88-8 | | | |

U.S. EPA Label Information

EPA Pesticide Registration Number Not regulated

16. OTHER INFORMATION

Issuing Date21-Apr-2009Revision Date08-Jan-2016

Revision Note

No information available

Disclaimer

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. VWR International and its Affiliates shall not be held liable for any damage resulting from handling.

End of Safety Data Sheet