



ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

AMRESCO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

Silver Nitrate

Code: 0377

Chemical Formula:	AgNO ₃	Manufacture Date:	(batch specific)
Molecular Weight:	169.87	Expiration/Reassay Date:	(batch specific)
CAS #:	7761-88-8		
Appearance:		Storage:	
Colorless or white solid.		Grade:	ACS GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Chloride	<= 0.0005 %	PASS
Clarity of Solution	PASS	PASS
Copper	<= 0.0002 %	PASS
Free Acid Test	PASS	PASS
Iron	<= 0.0002 %	PASS
Lead	<= 0.001 %	PASS
Purity	>= 99.0 %	PASS
Substances not Precipitated by HCl	<= 0.01 %	PASS
Sulfate	<= 0.002 %	PASS

Spec Set: 0377ACS

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____

Title:

Date Printed:

03/09/2016

Page 1 of 1