



ISO 9001 CERTIFIED
 ISO 13485 CERTIFIED

AMRESCO LLC
 28600 Fountain Parkway
 Solon, Ohio USA 44139
 440/349-1199 FAX: 440/349-1182
 www.amresco-inc.com
 Email: info@amresco-inc.com

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

Crystal Violet

Code: 0528

Chemical Formula:	C ₂₅ H ₃₀ N ₃ Cl	Manufacture Date:	(batch specific)
Molecular Weight:	408.00	Expiration/Reassay Date:	(batch specific)
CAS #:	548-62-9		
Appearance:		Storage:	
Green to dark green iridescent powder		Grade:	USP GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Arsenic	<= 0.001 %	PASS
Chromatographic Purity Test	PASS	PASS
Dye Content	96.0 - 100.5	PASS
Identification	PASS	PASS
Insolubles (Alcohol)	<= 1.0 %	PASS
Lead	<= 0.003 %	PASS
Residue on Ignition	<= 1.5 %	PASS
Water	<= 7.5 %	PASS
Zinc	<= 0.05 %	PASS

Spec Set: 0528USP

Analysis may have been rounded to significant digits in specification limits.
 Product meets analytical specifications of the grades listed.

Residual Solvents, if applicable will be reported for actual batch number only.

Internal ID #: 269

Signature: _____
 Title:

Date Printed: 05/04/2016
 Page 1 of 1