

# SAFETY DATA SHEET

Issuing Date 19-May-2009 Revision Date 08-Jan-2016 Revision Number 2

# 1. PRODUCT AND COMPANY IDENTIFICATION

**Product Identifier** 

Product Name 2-Mercaptoethanol

Other means of identification

Product Code(s) 0482 UN/ID no. 2966

Synonyms No information available

Recommended use of the chemical and restrictions on use

Recommended Use For Further Manufacturing Use Only.
Uses advised against Not for Human or Animal Drug Use

Details of the supplier of the safety data sheet

Company AddressManufacturer AddressDistributorVWR International, LLCVWR International, LLCAMRESCO INC.Radnor Corporate Center28600 Fountain Parkway6681 Cochran Road100 Matsonford RoadSolon, Ohio 44139SOLON, OHIO 44139

Radnor, PA 19087-8660

Company Phone Number 1.800.829.2805 E-mail Address info@amresco-inc.com

**Emergency Telephone Number** 

Emergency Telephone Number Chemtrec 1-800-424-9300

# 2. HAZARDS IDENTIFICATION

#### Classification

#### **OSHA Regulatory Status**

This chemical is considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200)

Acute Toxicity - Oral	Category 3
Acute Toxicity - Dermal	Category 2
Acute toxicity - Inhalation (Gases)	Category 3
Skin Corrosion/Irritation	Category 2
Serious Eye Damage/Eye Irritation	Category 1
Skin Sensitization	Category 1
Specific target organ toxicity (repeated exposure)	Category 2
Flammable liquids	Category 4

#### Label elements

### **Emergency Overview**

# Danger

#### Hazard statements

H301 - Toxic if swallowed

H310 - Fatal in contact with skin

H315 - Causes skin irritation

H317 - May cause an allergic skin reaction

H318 - Causes serious eye damage

H331 - Toxic if inhaled

H373 - May cause damage to organs through prolonged or repeated exposure

H227 - Combustible liquid



Appearance Clear Colorless

Physical State Liquid

**Odor** This material has a strong objectionable odor

# **Precautionary Statements - Prevention**

Wash face, hands and any exposed skin thoroughly after handling

Do not eat, drink or smoke when using this product

Do not get in eyes, on skin, or on clothing

Wear protective gloves/protective clothing/eye protection/face protection

Use only outdoors or in a well-ventilated area

Contaminated work clothing should not be allowed out of the workplace

Do not breathe dust/fume/gas/mist/vapors/spray

Keep away from heat/sparks/open flames/hot surfaces. - No smoking

Keep cool

#### **Precautionary Statements - Response**

Specific treatment (see .? on this label)

Specific measures (see .? on this label)

Specific treatment (see .? on this label)

Get medical advice/attention if you feel unwell

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing

Immediately call a POISON CENTER or doctor/physician

IF ON SKIN: Gently wash with plenty of soap and water

Immediately call a POISON CENTER or doctor/physician

Remove/Take off immediately all contaminated clothing

Wash contaminated clothing before reuse

If skin irritation or rash occurs: Get medical advice/attention

IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing

Call a POISON CENTER or doctor/physician

IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician

Rinse mouth

In case of fire: Use CO2, dry chemical, or foam for extinction

#### **Precautionary Statements - Storage**

Store locked up

Store in a well-ventilated place. Keep container tightly closed

# **Precautionary Statements - Disposal**

Dispose of contents/container to an approved waste disposal plant

#### Hazards not otherwise classified (HNOC)

Not regulated

#### Other Information

Very toxic to aquatic life with long lasting effects

#### 3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical Name	CAS-No	EC No.	Weight %	Trade Secret
2-Mercaptoethanol	60-24-2	200-464-6	95-100	Not applicable

# 4. FIRST AID MEASURES

First Aid Measures

**Eye Contact** Rinse thoroughly with plenty of water for at least 15 minutes, lifting lower and upper eyelids.

Consult a physician.

**Skin Contact**Wash off immediately with soap and plenty of water removing all contaminated clothes and

shoes.

**Inhalation** Move to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give

oxygen.

**Ingestion** Clean mouth with water and afterwards drink plenty of water. Never give anything by mouth

to an unconscious person.

Most important symptoms and effects, both acute and delayed

Most Important Symptoms/Effects No information available.

Indication of any immediate medical attention and special treatment needed

Notes to Physician Treat symptomatically.

#### 5. FIRE-FIGHTING MEASURES

### **Suitable Extinguishing Media**

Dry chemical, CO<sub>2</sub>, water spray or regular foam.

Unsuitable Extinguishing Media No information available.

#### **Specific Hazards Arising from the Chemical**

No information available.

**Hazardous Combustion** 

Carbon oxides. Sulfur oxides.

Products

**Explosion Data** 

Sensitivity to Mechanical Impact None. Sensitivity to Static Discharge None.

### **Protective Equipment and Precautions for Firefighters**

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

#### 6. ACCIDENTAL RELEASE MEASURES

### Personal precautions, protective equipment and emergency procedures

**Personal protection** Avoid contact with skin, eyes and clothing. Ensure adequate ventilation.

**Environmental Precautions** 

**Environmental Precautions** See Section 12 for additional Ecological information.

#### Methods and material for containment and cleaning up

**Methods for Containment** Prevent further leakage or spillage if safe to do so.

Methods for Cleaning Up Pick up and transfer to properly labeled containers. Ventilate area and wash spill site after

material pickup is complete. Cover with an activated carbon adsorbent.

# 7. HANDLING AND STORAGE

**Precautions for Safe Handling** 

Handling Handle in accordance with good industrial hygiene and safety practice.

Conditions for safe storage, including any incompatibilities

Storage Keep container tightly closed in a dry and well-ventilated place.

Incompatible Products Strong oxidizers. Metals.

# 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Control parameters

Exposure Guidelines This product does not contain any hazardous materials with occupational exposure limits

established by the region specific regulatory bodies.

Chemical Name	ACGIH TLV	OSHA PEL	NIOSH IDLH
2-Mercaptoethanol	-	-	-
60-24-2			

**Appropriate engineering controls** 

Engineering Measures Showers

Eyewash stations Ventilation systems.

Individual protection measures, such as personal protective equipment

**Eye/Face Protection** Safety glasses with side-shields.

**Skin and Body Protection** Wear protective gloves/clothing.

Respiratory Protection If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved

respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be

Remarks • Method

provided in accordance with current local regulations.

**Hygiene Measures** Handle in accordance with good industrial hygiene and safety practice.

#### 9. PHYSICAL AND CHEMICAL PROPERTIES

### Information on basic physical and chemical properties

Physical State Liquid

Appearance Clear Colorless Odor This material has a strong

objectionable odor

Color No information available Odor Threshold No information available

<u>Property</u> <u>Values</u>

**pH** No information available

Melting point/freezing point -100 deg °C Boiling Point/Range 157.2 °C Flash Point (High in °C) 73.9 °C

Evaporation Rate
No information available
No information available

Flammability Limit in Air
Upper flammability limit:

Upper flammability limit:
Lower flammability limit:
Vapor pressure
Vapor Density

No information available
No information available
No information available
No information available

Specific Gravity 1.14

### Product Code(s)0482 - 2-Mercaptoethanol

**Water Solubility** Very soluble in water Solubility in other solvents No information available Partition coefficient No information available 295 °C / 563.00 °F **Autoignition Temperature Decomposition Temperature** No information available Kinematic viscosity No information available Dynamic viscosity No information available **Explosive Properties** No information available **Oxidizing Properties** No information available

Other Information

Softening Point No information available

Molecular Weight 78.13

VOC Content No information available

Density 1.14

Bulk Density No information available

# 10. STABILITY AND REACTIVITY

#### Reactivity

No data available

# **Chemical Stability**

Stable under recommended storage conditions.

#### **Possibility of Hazardous Reactions**

None under normal processing.

**Conditions to Avoid** 

Extremes of temperature and direct sunlight.

**Incompatible Materials** 

Strong oxidizers. Metals.

# **Hazardous Decomposition Products**

Carbon oxides. Sulfur oxides.

# 11. TOXICOLOGICAL INFORMATION

# Information on likely routes of exposure

#### **Product Information**

**Inhalation** Toxic by inhalation.

**Eye Contact** Severely irritating to eyes.

**Skin Contact** May be fatal if absorbed through skin.

**Ingestion** Toxic if swallowed.

Chemical Name	LD50 Oral	LD50 Dermal	LC50 Inhalation
2-Mercaptoethanol	= 244 mg/kg (Rat)	112 - 224 mg/kg (Rabbit) = 150	-
60-24-2		μL/kg (Rabbit)	

#### Information on toxicological effects

**Symptoms** No information available.

# Delayed and immediate effects as well as chronic effects from short and long-term exposure

SensitizationNo information available.Mutagenic EffectsNo information available.CarcinogenicityNo information available.

# Product Code(s)0482 - 2-Mercaptoethanol

\_\_\_\_\_

Chemical Name	ACGIH	IARC	NTP	OSHA
2-Mercaptoethanol	-	-	-	=
60-24-2				

Reproductive Toxicity
STOT - single exposure
STOT - repeated exposure
Aspiration hazard

No information available.
No information available.
No information available.

Numerical measures of toxicity - Product Information

The following values are calculated based on chapter 3.1 of the GHS document .

# 12. ECOLOGICAL INFORMATION

#### **Ecotoxicity**

0% of the mixture consists of components(s) of unknown hazards to the aquatic environment

Chemical Name	Algae/aquatic plants	Toxicity to Fish	Crustacea
2-Mercaptoethanol	12: 72 h Desmodesmus subspicatus	46 - 100: 96 h Leuciscus idus mg/L	1.52: 48 h Daphnia magna mg/L
60-24-2	mg/L EC50	LC50 static	EC50

#### Persistence and Degradability

No information available.

#### **Bioaccumulation/Accumulation**

No information available.

Chemical Name	Partition coefficient
2-Mercaptoethanol	-0.056
60-24-2	

Other Adverse Effects No information available

# 13. DISPOSAL CONSIDERATIONS

#### Waste treatment methods

Waste Disposal Method Dispose of material in accordance with all federal, state, and local regulations.

**Contaminated Packaging** Do not re-use empty containers.

Chemical Name	RCRA - Halogenated Organic Compounds	RCRA - P Series Wastes	RCRA - F Series Wastes	RCRA - K Series Wastes
2-Mercaptoethanol 60-24-2	-	-	-	-

Chemical Name	California Hazardous Waste Status
2-Mercaptoethanol	-
60-24-2	

# 14. TRANSPORT INFORMATION

<u>DOT</u>

UN/ID no. 2966

Proper shipping name THIOGLYCOL

Hazard Class 6.1 Packing Group

IATA

**UN/ID** no. 2966

Proper shipping name THIOGLYCOL

\_\_\_\_\_

Hazard Class 6.1 Packing Group

# 15. REGULATORY INFORMATION

**International Inventories** 

**TSCA** Complies **DSL/NDSL** Complies **EINECS/ELINCS** Complies **ENCS** Complies Complies **IECSC** Complies **KECL PICCS** Complies Complies **AICS** 

#### Legend:

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS - European Inventory of Existing Commercial Chemical Substances/EU List of Notified Chemical Substances

**ENCS** - Japan Existing and New Chemical Substances

**IECSC** - China Inventory of Existing Chemical Substances

**KECL** - Korean Existing and Evaluated Chemical Substances

PICCS - Philippines Inventory of Chemicals and Chemical Substances

AICS - Australian Inventory of Chemical Substances

#### U.S. Federal Regulations

#### **SARA 313**

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

Component	SARA 313 - Threshold Values %
2-Mercaptoethanol	-
60-24-2 ( 95-100 )	

#### SARA 311/312 Hazard Categories

Acute Health Hazard
Yes
Chronic Health Hazard
Yes
Fire Hazard
Yes
Sudden Release of Pressure Hazard
No
Reactive Hazard
No

### **Clean Water Act**

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42)

Chemical Name	CWA - Reportable Quantities	CWA - Toxic Pollutants	CWA - Priority Pollutants	CWA - Hazardous Substances
2-Mercaptoethanol 60-24-2	-	-	-	-

#### CERCLA

This material, as supplied, does not contain any substances regulated as hazardous substances under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302) or the Superfund Amendments and Reauthorization Act (SARA) (40 CFR 355). There may be specific reporting requirements at the local, regional, or state level pertaining to releases of this material

Chemical Nam	e Hazardous S	Substances RQs	Extremely Hazardous Substances RQs	RQ
2-Mercaptoethar 60-24-2	nol	-	-	-

# **U.S. State Regulations**

# **California Proposition 65**

This product does not contain any Proposition 65 chemicals.

Chemical Name	California Prop. 65	
2-Mercaptoethanol - 60-24-2	-	

### **U.S. State Right-to-Know Regulations**

Chemical Name	New Jersey	Massachusetts	Pennsylvania
2-Mercaptoethanol	X	X	X
60-24-2			

#### **U.S. EPA Label Information**

EPA Pesticide Registration Number Not regulated

# **16. OTHER INFORMATION**

Issuing Date19-May-2009Revision Date08-Jan-2016

**Revision Note** 

No information available

**Disclaimer** 

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. VWR International and its Affiliates shall not be held liable for any damage resulting from handling.

**End of Safety Data Sheet**