

**AMRESCO LLC**

28600 Fountain Parkway  
 Solon, Ohio USA 44139  
 440/349-1199 FAX: 440/349-1182  
 www.amresco-inc.com  
 Email: info@amresco-inc.com

ISO 9001 CERTIFIED  
 ISO 13485 CERTIFIED

**CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS****Coomassie Brilliant Blue R-250****Code: 0472**

Chemical Formula:	C45H44N3O7S2Na	Manufacture Date:	(batch specific)
Molecular Weight:	825.99	Expiration/Reassay Date:	(batch specific)
CAS #:	6104-59-2		
Appearance:	Blue to red crystalline powder	Storage:	
		Grade:	ULTRA PURE GRADE
		Additional Product Description:	
		Cl#	42660

**Additional Information**

TEST	SPECIFICATION	DISPOSITION
Dye Content (Dry Basis)	>= 65 %	PASS
Em (550nm, Water, Dry Basis)	>= 33000	PASS
Identification	PASS	PASS
Loss on Drying	<= 5.0 %	PASS
Solubility (0.020g/L, Water)	PASS	PASS

Spec Set: 0472UP

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: \_\_\_\_\_

Title:

Date Printed:

03/08/2016

Page 1 of 1