

**AMRESCO LLC**

28600 Fountain Parkway
 Solon, Ohio USA 44139
 440/349-1199 FAX: 440/349-1182
 www.amresco-inc.com
 Email: info@amresco-inc.com

ISO 9001 CERTIFIED
 ISO 13485 CERTIFIED

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS**SODIUM PHOSPHATE, DIBASIC, HEPTAHYDRATE****Code: 0348**

Chemical Formula:	Na ₂ HPO ₄ ·7H ₂ O	Manufacture Date:	(batch specific)
Molecular Weight:	268.07	Expiration/Reassay Date:	(batch specific)
CAS #:	7782-85-6		
Appearance:		Storage:	Room Temperature
Fine white crystalline powder		Grade:	ACS GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Chloride	<= 0.001 %	PASS
Heavy Metals (as Pb)	<= 0.001 %	PASS
Insolubles	<= 0.005 %	PASS
Iron	<= 0.001 %	PASS
pH (5%, Water) @25C	8.7 - 9.3	PASS
Purity (as is)	98.0 - 102.0 %	PASS
Sulfate	<= 0.005 %	PASS

Spec Set: 0348ACS

Analysis may have been rounded to significant digits in specification limits.
 Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____

Title:

Date Printed:

03/08/2016

Page 1 of 1