



AMRESCO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

Separating Gel 4X Buffer

Code: M197

Chemical Formula:	Not Applicable	Manufacture Date:	(batch specific)
Molecular Weight:		Expiration/Reassay Date:	(batch specific)
CAS #:	Not Applicable		
Appearance:		Storage:	Room Temperature
Colorless, haze-free liquid.		Grade:	PROTEOMICS GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Color Test	PASS	PASS
Conductivity @25C	14.0 - 17.5 mmhos	PASS
pH @25C	8.7 - 8.9	PASS
Protease	NONE	PASS
Titration	1.4 - 1.6 M - Molar	PASS

Spec Set: M197PTM

Analysis may have been rounded to significant digits in specification limits.
Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____

Title:

Date Printed: 08/11/2016

Page 1 of 1