



**PCR 8-Strip Tube link with Cap**  
Certificate of Quality

**Cat. No. :** PT-G-0208-LN-C

**Product Description:** 0.2ml PCR 8-Strip Tube link with Caps

**Lot No. :** F1213-247094-0208LNC

**Material :** Virgin Polypropylene

**Quantity :** 60 strips/pk

**Expiry Date:** 2022/08/31

All Protech 0.2ml PCR 8-Strip Tube link with Cap (PT-G-0208-LN-C) are manufactured in a controlled area and are subject to regular quality inspections within statistical process control to meet non-detectable DNase, RNase, human DNA and Pyrogen.

The Protech Quality Management System is certified according to the international standards ISO 9001:2008.

Protech 0.2ml PCR 8-Strip Tube link with Cap (PT-G-0208-LN-C) are made of FDA and USP VI standard virgin Polypropylene of highest purity and quality. Our material suppliers do not use or intentionally incorporate the following agents into the materials

Protech is using for the production of Tubes:

- Slip agents (including oleamide, erucamide, stearamide)
- Biocides (including di(2-hydroxyethyl)methyldodecylammonium salts (DiHEMDA))
- Plasticizers (phthalates)

**Approved by:** Danny Hsu

**Title:** **QC SUPERVISOR**

**Date:** **2020.03.20**

For Research Using Only

Please do not hesitate to contact us if you have any questions.

Manufactured for and distributed by Protech Technology Enterprise Co., Ltd.

TEL: +886-2-2655-7677 / FAX: +886-2-2655-7601 / Toll Free: 0800-231-530

E-MAIL: service@bio-protech.com.tw; tech@bio-protech.com.tw

WEB: <http://www.bio-protech.com.tw>