

R.T.U X-Gal Solution

Cat No: PT-BL005

Size: 2mL

Concentration: 50 mg/ mL

Description:

Protech PT-BL005 X-Gal Solution is a ready-to-use, stabilized, non-toxic solution. The X-Gal is a non-inducing chromogenic substrate of β -galactosidase that hydrolyzes X-Gal, resulting in the formation of an intense blue precipitate. The X-Gal solution is most frequently used in combination with IPTG for blue/white colony screening to distinguish recombinants (white) from non-recombinants (blue).

Storage buffer:

The RTU X-Gal solution is supplied in a non-hazardous and non-toxic solvent.

Storage:

Stable for a year when stored at -20°C.

*****Research Use Only*****

Please do not hesitate to contact us if you have any questions.

Manufactured for and distributed by Protech Technology Enterprise Co., Ltd.

TEL : +886-2-2655-7677 / FAX : +886-2-2655-7601 / Toll Free : 0800-231-530

E-MAIL: service@bio-protech.com.tw ; tech@bio-protech.com.tw

WEB: <http://www.bio-protech.com.tw>