

UltraPure Water
Certificate of Analysis

Cat. No. : PT-P570

Size: 1 L

PT-P570-50

50 mL

Lot. No. : JD292404

Store: Room Temperature

Re-assay Date: 2024.4.30

Product Description :

Nuclease Free Water is a high quality reagent grade water produced by a professional water purification system. After production the water is sterile filtered and dispensed into sterile RNase free containers. The water is tested for DNase, RNase and proteinase contamination through the standard method.

The Nuclease Free Water is indicated for reagent preparation in molecular biology, nucleic acid research, diagnostics, tissue culture, in vitro fertilization, residue analysis, HPLC, spectrophotometry, and for many other application requiring exceptional quality water.

Quality Control :

DNase Activity ----- Non Detected

RNase Activity ----- Non Detected

Proteinase Activity ----- Non Detected

Sterility ----- Pass

Approved by: Jiang, Chao-Sheng

Title: **QC SUPERVISOR**

Date: **2022.07.21**

For Research Using Only

Please do not hesitate to contact us if you have any questions.

Manufactured for and distributed by Protech Technology Enterprise Co., Ltd.

TEL: +886-2-2655-7677 / FAX: +886-2-2655-7601 / Toll Free: 0800-231-530

E-MAIL: service@bio-protech.com.tw; tech@bio-protech.com.tw

WEB: <http://www.bio-protech.com.tw>