

UltraPure DEPC Water

Certificate of Analysis

Catalog No. : PT-P560-500

Lot No. : IL282312

Quantity : 500mL

Storage: Room Temperature

Re-assay Date: 2023.12.31

Product Description:

UltraPure DEPC Water is high quality reagent grade water produced by a professional water purification system. The water is treated with diethylpyrocarbonate (DEPC) and autoclaved to remove the residual DEPC. After production the UltraPure DEPC Water is dispensed into sterile RNase free containers and ready for use. It is tested for DNase, RNase and proteinase contamination through the standard method.

The UltraPure DEPC Water can be used for reagent preparation in molecular biology, nucleic acid research, diagnostics, and many other applications that require RNase- free water.

Quality Control :

DNase Activity	Not Detected
RNase Activity	Not Detected
Proteinase Activity	Not Detected
Sterility	Pass

Approved by: Jiang, Chao-Sheng

Title: **QC SUPERVISOR**

Date: **2022.04.27**

For Research Using Only

Please do not hesitate to contact us if you have any questions.

Manufactured for and distributed by Protech Technology Enterprise Co., Ltd.

TEL: +886-2-2655-7677 / FAX: +886-2-2655-7601 / Toll Free: 0800-231-530

E-MAIL: service@bio-protech.com.tw; tech@bio-protech.com.tw

WEB: <http://www.bio-protech.com.tw>