



ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

AMRESKO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

TRIS 1.5M PH8.8

Code: J831

Chemical Formula:	Not Applicable	Manufacture Date:	(batch specific)
Molecular Weight:		Expiration/Reassay Date:	(batch specific)
CAS #:	Not Applicable		
Appearance:		Storage:	Room Temperature
Haze-free, colorless liquid		Grade:	BIOTECHNOLOGY GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Conductivity @25C		PASS
pH @25C	8.7 - 8.9	PASS
Purity	16.96 - 19.38 %	PASS

Spec Set: J831BIO

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____

Title:

Date Printed:

06/23/2016

Page 1 of 1