

2020

PDM2020011
即日起至2020年底

波仕特年終清倉特賣會

滿\$5,000 四送一

RNase/DNase/EtBr
Terminator (R475-100mL) X1
Timer X1
7-11商品卡 NTD\$100 X2
威秀電影票 X1

滿\$10,000 三送一

7-11商品卡 NTD\$100 X4
NaCl (0241-1KG) X2
小米手環 X1

滿\$20,000
20% OFF!

核酸電泳

Catalog No.	產品名稱 / 包裝	特性
VWR-0478-2PK	10X TBE Ready-Pack, 2pk/bag	DNA 電泳用 TBE 緩衝液
VWR-0912-2PK	25X TAE Ready-Pack, 2pk/bag	DNA 電泳用 TAE 緩衝液
VWR-X328-10ML	Ethidium Bromide, 10mg/ml, 10ML	DNA 染色劑, 毒化物
VWR-N313-1ML/KIT	EZ-Vision Three DNA Loading Buffer, 1or 5ml	DNA 安全染劑, 含三種可見 Dyes
VWR-N650-KIT	EZ-Vision Two DNA Loading Buffer, 5ml	DNA 安全染劑, 含二種可見 Dyes
VWR-N472-KIT	EZ-Vision One DNA Loading Buffer, 5ml	DNA 安全染劑, 含一種可見 Dyes
VWR-E526-500ML	MOPS 10X Liquid, 500ml	RNA 電泳用緩衝液

細菌培養/抗生素

VWR-J859-500G	Tryptone, 500g	胰水解蛋白
VWR-J869-500G	Terrific Broth, 500g	細菌培養液, 須外加Glycerol
VWR-0339-25G	Ampicillin Sodium Salt, 25g	Working conc. 50~100ug/ml
VWR-0408-10G	Kanamycin Sulfate, 10g	Working conc. 25~50ug/ml
VWR-J593-25MG	Puromycin, Dihydrochloride, 25mg	Working conc. 1~10ug/ml
VWR-E859-5G	G418 Sulfate, 5g	Working conc. 100~200ug/ml

核酸/蛋白質純化

VWR-N653-100ML	RIPA Lysis Buffer, 100ml	蛋白質裂解液
VWR-0706-100MG	Proteinase K, 100mg	蛋白酶 K, >30U/mg
VWR-M250-1ML	Protease Inhibitor Cocktail, Mammalian, 1ml	哺乳細胞用抗蛋白酶組合液
VWR-0527-100G	Imidazole, 100g	純化 his-tag 蛋白質

蛋白質電泳

Catalog No.	產品名稱 / 包裝	特性
VWR-E344-500ML	Acryl/Bis 29:1 30%, 500ml	蛋白質/Native DNA/RNA電泳
VWR-E347-500ML	Acryl/Bis 37.5:1 30%, 500ml	蛋白質電泳, SDS-PAGE 用
VWR-0486-25G/100G	Ammonium Persulfate (APS), 25g/100g	10% APS solution, SDS-PAGE 用
VWR-0761-25ML	TEMED, 25ml	SDS-PAGE 搭配APS使用
VWR-M337-5ML	Laemmli Loading Buffer ,4X, 5ml	蛋白質電泳 Loading Dye
VWR-0837-500ML	SDS Solution 20%, 500ml	SDS-PAGE 用
VWR-0227-250G/500G/1KG	SDS, 250g/500g/1kg	SDS-PAGE 用

西方轉染



VWR-E116-100ML	BCIP/NBT Solution, 100ml	Western blot 呈色試劑
VWR-N552-1L	Gentle Stripping Buffer, 1L	溫和型抗體去除緩衝液

有機溶劑

VWR-0231-500ML	DMSO, 500ml	ACS grade
VWR-N182-5X10ML	DMSO, Ultrapure, 5x10ml	Cell culture 等級
VWR-0323-4L	METHANOL ANHYDROUS, 4L	>99.8% Ultra pure grade
VWR-0918-1L/-4L	Isopropyl Alcohol, 1L /4L	異丙醇; propan-2-ol
VWR-0883-400ML	Phenol:Chloroform, pH6.7/8.0, 400ml	PCI (25:24:1) 可用於蛋白質純化
VWR-0945-400ML	Phenol, Saturated, pH6.6/7.9, 400ml	可調至 pH 7.9, 用於DNA純化
VWR-0981-400ML	Phenol, Saturated, pH4.5, 400ml	用於RNA純化

緩衝液



VWR-J373-4L	PBS 10X, pH:7.4, 4L	PBS 緩衝液
VWR-E404-100/-200TABS	PBS Tablets, 100ml/tablets, 100/200 tabs	PBS 緩衝液, 錠狀
VWR-0780-2PK/-50L	PBS, Ready-Pack, 2pk/bag or 50L	PBS 緩衝液, cell culture 等級
VWR-0788-2PK	TBS Buffer, Ready-Pack, 2pk/bag	TBS 緩衝液, 1L 20X buffer/pk

常用化學品



VWR-0188-1KG	D-Glucose, Anhydrous, 1kg	VWR-0556-500G	CaCl ₂ , Dihydrate, 500g
VWR-0281-5G	DTT, 5g 還原劑	VWR-0583-500G	NaOH, Pellets, 500g
VWR-0288-100G	Magnesium Chloride, 100g	VWR-0781-500G	KH ₂ PO ₄ , Anhydrous, 500g
VWR-0332-100G	Albumin Bovine(BSA), 100g	VWR-0793-1G	MTT, 1g 細胞活性呈色劑
VWR-0335-1KG/5KG	Sucrose, 1kg/5kg	VWR-0865-1KG	Sodium Bicarbonate, 1kg
VWR-0404-500G	Na ₂ HPO ₄ , Anhydrous, 500g	VWR-E174-25G	DEPC, 25g(20.8-22.7ML)

