

NACALAI TESQUE, INC.

Nijo Karasuma, Nakagyo-ku, Kyoto 604-0855, JAPAN

Product No : 11585 Grade : SP for Cell Culture

Product Name : DMEM (4.5g/l Glucose) with HEPES, without L-Gln and Sodium Pyruvate, liquid

MF :

FW :

Test item	Unit	Specified value	Test method
pH (25deg.C)		7.0 – 7.5	JIS K8001
Osmolality	mOsm/kg	295 – 335	Cryoscopic method
Endotoxin	EU/mL	max. 0.2	Chromogenic synthetic substrate method
Mycoplasma		passes test	MycoAlert assay
Sterility test		passes test	Culture method