



October 01 2025

NACALAI TESQUE, INC.

Kyoto Factory

17 Ishibashi, Kaide-cho,
Muko, Kyoto 617-0004
TEL : +81- (0) 75-251-1730

Product number Grade Lot number Package size
13408-64 SP L5M1005 5X5ML

Product name
Dimethyl Sulfoxide

Test Item	Unit	Specification	Result
Purity	%	min. 99.5	100.0
Density (20deg.C)	g/mL	1.099 - 1.103	1.101
Refractive index n ₂₀ /D		1.477 - 1.480	1.479
Water	%	max. 0.2	max. 0.2
Endotoxin	EU/g	max. 0.5	max. 0.5
Sterility test		passes test	passes test
Suitability test		passes test	passes test

Remarks

k. Goya

SIGNATURE