

Date: 2026.02

Version 26_01

Fungi Clean (200X)

Cat No : PT-FC200-500ML

Size : 500mL

Description

Fungi Clean is a specialized formula designed for water treatment. It exhibits broad-spectrum antimicrobial activity against Gram-positive bacteria, Gram-negative bacteria, yeast, mycoplasma, and algae.

Specifically, Fungi Clean is non-toxic to higher eukaryotic cells and is therefore more environmentally friendly than similar products when used in water reservoirs inside incubators.

In addition to its antimicrobial properties, Fungi Clean can also prevent limescale and rust formation when used in water baths.

Protocol

Dilute 100-200 folds with distilled water and add diluted fluid into incubator water reservoirs or water bath.

Renew fluid every 4 weeks, or the Fungi Clean fluid blue color is decayed.

- Fungi Clean is stored at room temperature for 2 years.
- For research use only. NOT FOR USE IN FOOD, DRUG OR MEDICAL DEVICE MANUFACTURING.

Research Use Only

Please do not hesitate to contact us while you have any questions.

Manufactured for and distributed by Protech Technology Enterprise Co., Ltd.

TEL : +886-2-2655-7677 / FAX : +886-2-2655-7601 / Toll Free : 0800-231-530

E-MAIL: service@bio-protech.com.tw ; tech@bio-protech.com.tw

WEB: <http://www.bio-protech.com.tw>