



August 27 2026

NACALAI TESQUE, INC.

Kyoto Factory

17 Ishibashi, Kaide-cho,
Muko, Kyoto 617-0004
TEL : +81- (0) 75-251-1730

Product number Grade Lot number Package size
12679-54 SP V6M9713 100ML

Product name
Accutase (TM)

Expiration date Jul 31 2028

Test Item	Unit	Specification	Result
Appearance		Clear, light pink solution	passes test
pH		7.2 - 7.5	7.3
Trypsin activity	EAU/min/mL	500 - 720	667
Chymotrypsin/Elastase activity	EAU/min/mL	500 - 900	564
Collagenase Type I activity	RFU/min/mL	350 - 750	414
Cell detachment activity		passes test	passes test
Endotoxin	EU/mL	max. 1	max. 1
Sterility test		passes test	passes test

Remarks

M. Kawashimo

SIGNATURE