ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

AMRESCO LLC
28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

Sodium Bicarbonate

Code: **0865**

Chemical Formula:	$NaHCO_3$	Manufacture Date:	(batch specific)
Molecular Weight:	84.01	Expiration/Reassay Date:	(batch specific)
CAS #:	144-55-8		
Appearance:		Storage:	
Fine white crystalline powder		Grade:	ACS GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Ammonium	<= 0.0005 %	PASS
Calcium	<= 0.02 %	PASS
Chloride	<= 0.003 %	PASS
Heavy Metals (as Pb)	<= 0.0005 %	PASS
Insolubles	<= 0.015 %	PASS
Iron	<= 0.001 %	PASS
Magnesium	<= 0.005 %	PASS
Phosphate	<= 0.001 %	PASS
Potassium	<= 0.005 %	PASS
Purity	99.7 - 100.3 %	PASS
Sulfur Compounds	<= 0.003 %	PASS

Spec Set: 0865ACS

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____ **Date Printed:** 06/08/2016

Title: